

Electrochemical Impedance Spectroscopy (EIS) for Advanced Battery Diagnostics

> Prepared by: Sujan Shrestha, PhD & Mark Sholin



Electrochemical Impedance Spectroscopy

Electrochemistry

AC Resistance

 $\mathbf{Z} = \frac{V_{AC}}{i_{AC}}$

IN OTHER WORDS,

EIS is a **non-invasive**, typically **rapid** measurement technique that reveals a wealth of information about conditions inside a battery



Range of Frequency

Three Most Common EIS Applications for Batteries

Battery Testing
 Battery Management
 Battery R&D

Raw Information

- Internal Resistance
- Charge-transfer Resistance
- Mass-transfer Resistance
- Interfacial Resistance
- State of Charge (SOC)
- State of Health (SOH)

Useful Applications

- Battery Characterization
- Battery Management System
- New Battery Development
- Failure Diagnosis



Steps to run a standard EIS experiment

- Find an EIS-capable instrument, typically a potentiostat, and connect it to the battery being tested
- 2. Decide to use Potentiostatic EIS or Galvanostatic EIS on the battery. *Pot. = Voltage input, Galv. = Current input*
- **3.** Configure input AC signal the potentiostat will apply to the battery (frequency range, amplitude, DC bias, etc.)



The potentiostat applies the AC input signal and records the resulting AC output signal for internal processing



Internal processing calculates: Frequency, Impedance Modulus (real & i), Phase, DC Current, and DC Voltage

Analyze the resulting data, typically in the form of a Nyquist Plot fitted to an equivalent circuit model





EIS Data is Typically Graphed Two Ways



Nyquist Plots are more often used because several Key Performance Indicators (KPIs) can be quickly derived through quick visual inspection

Ai

EIS Profile of a Typical Li-ion Battery

The **Nyquist Plot** of a typical Li-ion battery shows **four** prominent regions



Real Impedance (Ohm)



Frequency-dependent regions are correlated to many phenomena inside a battery

Very High Frequency Region Li-ion conductivity in electrolyte, electronic conductivity

High Frequency Region Li-ion transport through solidelectrolyte interphase

Mid Frequency Region Charge transfer at the electrodeelectrolyte interface

Low Frequency Region Li-ion diffusion in electrode and electrolyte



Equivalent Circuit Modeling (ECM) for EIS Data

Generic Electrochemical System



The behavior of electrochemical systems can be parameterized using the "language" of circuit-like elements (resistors, capacitors, inductors, etc.)



By fitting EIS data to an ECM, the electrochemical phenomena can be *quantified*

Examples:

R_s : Series resistance

R_{SEI} : Solid-electrolyte interface resistance

R_{cT}: Charge-transfer resistance

R_D: Diffusion resistance

Imaginary Impedance (Ohm) **High Frequency Region** Mid Frequency Region Very High Frequency Region Series Resistance SEI Resistance Charge-transfer Resistance **Diffusion Resistance** Rs R_{SEI} R_{CT} R_{D} Equivalent Circuit Model

State-of-Charge (SOC) and State-of-Health (SOH) values can be derived from the constants calculated by the equivalent circuit model



Calculating ECM Constants Needs Special Software



ECM software is available from most potentiostat vendors to perform analysis



Circuit Elements	Fitted Values
R1	152 m0hm
C1	5.53 F
R2	22.5 m0hm
C2	21.1 F
R3	67.6 m0hm
C3	27.3 F
R4	8.85 m0hm
W1	3.04 mDW

Ai

Series resistance (R_S) increases with cycle life, as shown by the larger x-axis intercept at higher cycles



Figure from: Jiang, et al., Energy Procedia 105 (2017)844-849



Charge-transfer resistance (R_{CT}) increases with SOC, as shown by taller & longer "humps" in the middle



Figure from: M. Li, Renewable Energy 100 (2017) 44-52



Pros & Cons of EIS Versus Traditional DC Techniques Pros Cons

- EIS is a non-invasive method (won't prematurely age cells, unlike some DC methods)
- EIS generates results in seconds or minutes, versus hours or days for charge-discharge tests
- EIS does not require as high-power of test
 instruments compared to DC-IR tests
- ECM reference data can be a "one-stop-shop" for KPIs like SOH, SOC, rather than running separate DC tests to obtain one KPI at a time

- EIS data analysis is generally more complicated than interpreting DC data
- ECMs can be inaccurate if reference data is not collected or analyzed properly
- Must carefully review EIS equipment specs to ensure they are compatible with your battery characteristics
- EIS is more sensitive to external sources of measurement noise

The ultimate solution is to have equipment capable of running both EIS and DC techniques. Most potentiostats (and some battery cyclers) can.



How To Incorporate **EIS** Into Your Testing Protocols

Identify the Key Performance Indicators (KPIs) that you want to measure (SOH, SOC, etc.) Run EIS on reference batteries operating with known KPIs in well-controlled conditions Develop Equivalent Circuit Models (ECMs) using EIS data from reference batteries

Measure KPIs of unknown batteries by running EIS on them & mapping data to ECMs

The accuracy and sensitivity of these KPI measurements depends upon the quantity and quality of reference data used to make the ECMs!





Ready to run EIS on your batteries but don't have the right tools?

Squidstat M NEXT GENERATION Electrochemical Testing Tools

Why Choose Us?

Research-grade precision at widely accessible prices

Powerful, compact hardware adds versatility & saves space

Increases workflow and saves times with intuitive software



<u>Click Here</u> to register for a FREE 30-Day Trial (USA only)

-@



FORWARD-THINKING HARDWARE. INTUITIVE SOFTWARE.

We're Devoted to Making Your Job EASIER



Mark Sholin Co-founder & President mark@admiralinstruments.com Office: +1 480 256 8706



Making the next generation of electrochemistry instruments **truly accessible worldwide**



Confidential